



Company:

Customer:

Date:

Basic price: **Eltego**
Max. capacity 120 kg

€ 3.558,-

▼ Please mark!

Seat	€
Seat width adjustable 400 - 540 mm	●
Seat height 485 mm	●
Seat depth 440 mm	●
Seat cushion	●
Standard seat black	●
Back	
Back support height 540 mm	●
Back support angle adjustable from 90° to 120°, mechanically	●
Back support adjustable by strap system	●
Head support incl. fixation	144,-
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Footplate angle adjustable	●
Motor and Battery	
Max. range 25 km	●
Speed 6 km/h	●
Batteries 2x 12V / 50Ah	●
Battery charger 24V / 4Ah	●

▼ Please mark!

Control device	€
swing-away control device, adjustable in length	●
Control device mounted left right	○
Tyres	
Rear wheel suspension	●
Rear wheel 12,5" x 2 1/4", PU	●
Front wheel 200 x 50 mm, PU	●
Mounting PU to pneumatic tires	79,-
Accessories	
Front and rear lights	●
Safety belt	●
Crutch holder	56,-
APINO Shopper	18,20
Tool set	●
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge