



ORDER FORM

TOGO
Power pushing aid

Company:

Customer:

Date:

Basic price: **Togo**

€ 2.180,-

Max. total capacity 180 kg

(including weight of pushing aid and wheelchair)

▼ Please mark!

Technical Data	€
Togo for seat width 37 - 52 cm	●
Togo for seat width 52 - 56 cm	201,80
Speed forward 1,0 up to 5,5 km/h 6 levels	●
Speed backwards approx. 1,5 km/h	●
Max. range approx. 12 km	●
Upward slope / downward slope max. 10%	●
Battery 24 V 11,5 Ah Standard	●
Operating voltage: 24 VDC	●
Battery charger	●
Drive wheels 200 x 50 solid rubber tires	●
Wheel lock electromagnetic	●
Software kit	89,60

▼ Please mark!

Accessories	€
Battery Power: 24 V 14,4 Ah	873,30
Adaptation for wheelchair mounting only with tool holder	221,40
Tool holder for wheelchair	45,20
Tool holder for multifunctional wheelchair Protego	121,60
Mounting aid for adaptation	189,50
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge