



ORDER FORM

CENTURO MINI 4-Wheel Scooter

Company:

Basic price: **Centuro Mini**

€ 1.504,20

Customer:

Max. capacity 136 kg

Date:

▼ Please mark !

Technical Data	€
Motor 6 km/h	●
- Range approx. 12 km	●
Battery 2x12 V / 12 Ah	●
Battery charger 3 A	●
Tyres 8" PU	●
Ground clearance 55 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21 %	●
Reversing width 1500 mm	●
Standard equipment	
Sprung suspension front and rear	●
Brakes electromagnetic	●
Handlebars ergonomic	●
Steering column adjustable	●
USB-Port	●
Cockpit with LCD display	●
Horn	●
Bicycle headlight	●

▼ Please mark !

Seat and back support	€
Swing-away seat, height-adjustable	
- Back support foldable	●
- Arm support angle adjustable	
Accessories	
Shopping basket removable	●
Doormat	●
Rear mirror left*	25,21
Rear mirror holder left	2,50
Rear mirror right*	25,21
Rear mirror holder right	2,50
Cane holder	92,44
Frame color	
White Anthracite	○

*just orderable with rear mirror holder

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge