



Company:

Customer:

Date:

▼ Please mark !

Model	€
Centuro S2 6 km/h - Range approx. 35 km - Max. capacity 160 kg	2.428,57
Centuro S2 10 km/h - Range approx. 30 km - Max. capacity 160 kg	2.932,77
Technical Data	
Battery 2x12 V / 45 Ah	●
Battery charger 5 A	●
Tyres 10" pneumatic	●
Ground clearance 80 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21 %	●
Reversing width 1950 mm	●
Standard equipment	
Sprung suspension front and rear	●
Electromagnetic brakes	●
Mechanical wheel-lock	●
Ergonomic handlebar	●
Adjustable steering column	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Rear mirror mounted left	●
Lighting according to German road traffic regulations	●
Authorized for public transportation	

▼ Please mark !

Seat and back support	€
Swing-away captain's seat, height-adjustable	
- Adjustable seat depth	●
- Foldable back support	
- Arm support, angle-adjustable	
Accessories	
Shopping basket removable	●
Doormat	●
Cup holder	●
Rear mirror mounted right*	25,21
Rear mirror holder right	2,50
License plate holder	24,10
Frame color	
White Anthracite	○

*just orderable with rear mirror holder

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge