



**Company:**

**Customer:**

**Date:**

▼ Please mark !

Model	€
Centuro S4   6 km/h - Range approx. 35 km - Max. capacity 160 kg	2.983,19
Centuro S4   10 km/h - Range approx. 30 km - Max. capacity 160 kg	3.605,04
<b>Technical Data</b>	
Battery 2x12 V / 55 Ah	<input type="radio"/>
Battery 2x12 V / 80 Ah	252,10*
Battery charger 5 A	<input checked="" type="radio"/>
Tyres 12"   pneumatic	<input checked="" type="radio"/>
Ground clearance 80 mm	<input checked="" type="radio"/>
Obstacle height 50 mm	<input checked="" type="radio"/>
Climbing ability 12°   21 %	<input checked="" type="radio"/>
Reversing width 2300 mm	<input checked="" type="radio"/>
<b>Standard equipment</b>	
Sprung suspension   front and rear	<input checked="" type="radio"/>
Electromagnetic brakes	<input checked="" type="radio"/>
Mechanical wheel-lock	<input checked="" type="radio"/>
Ergonomic handlebar	<input checked="" type="radio"/>
Adjustable steering column	<input checked="" type="radio"/>
USB-Port	<input checked="" type="radio"/>
Cockpit with LED display	<input checked="" type="radio"/>
Horn	<input checked="" type="radio"/>
Rear mirror   mounted left	<input checked="" type="radio"/>
Lighting according to German road traffic regulations	<input checked="" type="radio"/>
Electronic speed reduction	<input checked="" type="radio"/>

▼ Please mark !

Seat and back support	€
Swing-away captain's seat, height-adjustable - Adjustable seat depth - Foldable back support - Arm support, angle-adjustable	<input type="radio"/>
Swing-away captain's seat XL, height-adjustable - Adjustable seat depth - Foldable back support - Arm support, angle-adjustable	150,-*
<b>Accessories</b>	
Shopping basket   removable	<input checked="" type="radio"/>
Doormat	<input checked="" type="radio"/>
Cup holder	<input checked="" type="radio"/>
Rear mirror   mounted right**	25,21
Rear mirror holder   monted right	2,50
License plate holder	24,10
Cane holder	92,44
Transport box with lock   incl. fixture	250,-
Holder for standard rollators	216,-
<b>Frame color</b>	
White      Anthracite	<input type="radio"/>

\*Surcharge. Option price less price of standard equipment.  
This price differs from the spare part price.

\*\*just orderable with rear mirror holder

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment       Optional, no extra charge