



Company:

Basic price: **Centuro S7** | 15 km/h
Max. capacity 160 kg

€ 3.941,18

Customer:

Date:

▼ Please mark !

Technical Data	€
Motor 15 km /h - Reichweite ca. 45 km	●
Battery 2x12 V / 80 Ah	●
Battery charger 8 A	●
Tyres 13" pneumatic	●
Ground clearance 130 mm	●
Obstacle height 100 mm	●
Climbing ability 12° 21 %	●
Reversing width 2300 mm	●
Standard equipment	
Sprung suspension front and rear	●
Electromagnetic brakes	●
Mechanical wheel-lock	●
Ergonomic handlebar	●
Adjustable steering column	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Rear mirror mounted left and right	●
Lighting according to German road traffic regulations	●
Electronic speed reduction	●

▼ Please mark !

Seat and back support	€
Swing-away captain's seat, height-adjustable - Adjustable seat depth - Foldable back support - Arm support, angle-adjustable	○
Swing-away captain's seat XL, height-adjustable - Adjustable seat depth - Foldable back support - Arm support, angle-adjustable	150,-*
Accessories	
Shopping basket removable	●
Doormat	●
Cup holder	●
Licence plate holder	24,10
Cane holder	92,44
Transport box with lock incl. fixture	250,-
Holder for standard rollators	216,-
Frame color	
White Anthracite	○

*Surcharge. Option price less price of standard equipment.
This price differs from the spare part price.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge