



ORDER FORM

TRITON Multifunctional wheelchair

Company:

Basic price: **Triton without drum brake**

€ 1.999,-

Basic price: **Triton with drum brake**

€ 2.209,-

Customer:

Max. capacity 125 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm	
Variation I: 360 - 430	<input type="radio"/>
Variation II: 430 - 500	
Seat height 415 mm with seat cushion 520 mm	<input checked="" type="radio"/>
Seat depth 420 - 480 mm	<input checked="" type="radio"/>
Seat angle adjustable 0° - 15°	<input checked="" type="radio"/>
Medial support	92,70
Seat upholstery Dartex night grey	<input checked="" type="radio"/>
Back	
Comfort back support, shaped ergonomically	<input checked="" type="radio"/>
Back support height adjustable	<input checked="" type="radio"/>
Back support reclining from -3° up to +43°	<input checked="" type="radio"/>
Lateral trunk support swing away	<input checked="" type="radio"/>
Handgrip elevating, removable 850 - 1020 mm	<input checked="" type="radio"/>
Back support upholstery Dartex night grey	<input checked="" type="radio"/>
Head support adjustable in 3 directions	<input checked="" type="radio"/>
Head support extension	66,-
Sidepanels	
Sidepanels adjustable in height and depth, removable with upholstered arm support	<input checked="" type="radio"/>
Hemiplegic arm support	left right 298,70

▼ Please mark!

Leg support	€
Leg support elevating horizontally, with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" PU	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Accessories	
Anti-tip device Pair	87,50
Safety belt	26,80
Therapy tray	215,-
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge