



ORDER FORM

PROTEGO Multifunctional wheelchair

Company:

Basic price: **Protego without drum brake** € 1.990,-

Basic price: **Protego with drum brake** € 2.200,-

Customer:

Max. capacity 130 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 390 - 415 440 - 465 490 - 515	<input type="radio"/>
Seat height front 430 mm without cushion	<input checked="" type="radio"/>
Seat depth adjustable from 420 - 480 mm	<input checked="" type="radio"/>
Seat angle adjustable 0 - 15°	<input checked="" type="radio"/>
Medial support	92,70
Seat upholstery Dartex, night grey	<input checked="" type="radio"/>
Back	
Back support ergonomic shape with lateral support	<input checked="" type="radio"/>
Back support height 580 - 640 mm	<input checked="" type="radio"/>
Back support angle 0 - 44°	<input checked="" type="radio"/>
Trunk support swing away left right	94,80
Tool holder lateral trunk support	13,90
Back support adjustable by strap system	<input checked="" type="radio"/>
Back support upholstery Dartex, night grey	<input checked="" type="radio"/>
Handgrip elevating, removable	<input checked="" type="radio"/>
Sidepanels	
Sidepanels adjustable in height and depth, removable with arm support upholsteries Dartex upholstery	<input checked="" type="radio"/>
Hemiplegic arm support left right	298,70
Leg support	
Leg support elevating and angle adjustable with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>
Footplate only SW 44 cm	194,70

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" x 45 mm PU	<input checked="" type="radio"/>
Transit version, rear wheels 12,5" PU	210,-
Spoke protection, transparent 24"	145,-
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension	19,90
Accessories	
Togo power push aid incl. mounting set	2.523,-
Mounting aid for adaptation	189,50
Tilting aid	<input checked="" type="radio"/>
Safety belt	26,80
Therapy tray	215,-
Additional package Therapy tray with lateral trunk support	316,20
Oxygen/infusion holder 110 mm 160 mm	169,-
APINO Citybag	60,42
APINO Shopper	18,20
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge