



ORDER FORM

PYRO START PLUS SL

Light weight wheelchair

Company:

Customer:

Date:

Basic price: **Pyro Start Plus SL without drum brake** € 1.392,52
 Basic price: **Pyro Start Plus SL with drum brake** € 1.641,52
 Max. capacity 150 kg

▼ Please mark!

| Seat | | € |
|---|--|--------|
| Seat width in mm 430 460 490 520 | | ○ |
| Seat height front in mm 500 adjustable 525 - 600 mm | | ● |
| Seat depth 460 mm adjustable from 460 - 520 mm | | ● |
| Seat cushion Seat width 430-520, Seat depth 460-480 mm | | 30,- |
| Seat cushion Seat width 430-520, Seat depth 500-520 mm | | 34,- |
| Seat upholstery, Nylon | | ● |
| Back support | | |
| Back support height 460 mm adjustable 460 - 500 mm | | ● |
| Back support with lumbar curve | | ● |
| Back support adjustable by strap-system | | ● |
| Head support* incl. fixation | | 188,- |
| Handgrip elevating*, removable Pair | | 152,40 |
| Back support upholstery, Nylon | | ● |
| Sidepanels | | |
| Sidepanel removable with arm support, swing-away, elevating | | ● |
| Sidepanel elevating left right | | 108,20 |
| Sidepanel elevating with foldable arm pads left right | | 123,50 |
| Leg support | | |
| Swing-away leg support with individual footplates plastic, removable, adjustable in height | | ● |
| Leg support elongation 100 mm Pair | | 60,80 |
| Lower leg support elevating left right | | 107,60 |
| Amputee support left right | | 157,- |
| Footplates with heel straps | | ● |

* Head support and handgrips elevating cannot be mounted together

▼ Please mark!

| Tyres | | € |
|--|--|----------|
| Rear wheels 24" x 1 3/8" PU incl. quick-release axle | | ● |
| Front wheels 8" x 1 1/4", PU | | ● |
| Wheel base extension | | ● |
| Adapter block | | ● |
| Transit version, rear wheels 12,5" PU | | 210,- |
| Front wheels 200 x 50 mm. PU Pair | | 60,80 |
| Pneumatic tyres | | 79,- |
| Spoke protection, transparent 24" | | 145,- |
| One arm drive only without drum brake left right | | 325,- |
| Hand rim rubber-coated 24" Pair | | 110,- |
| Wheel locks | | |
| Knee lever wheel lock push to lock | | ● |
| Wheel lock lever extension | | 23,- |
| Knee lever wheel lock for one-arm drive left right | | 88,60 |
| Accessories | | |
| Togo power push aid incl. mounting set | | 2.446,60 |
| Mounting aid for adaptation | | 189,50 |
| Anti-tip device Pair | | 87,50 |
| Tilting aid | | ● |
| Safety belt | | 26,80 |
| Crutch holder | | 23,70 |
| Therapy tray | | 215,- |
| Oxygen/infusion holder and bracket | | 169,- |
| APINO Citybag | | 60,42 |
| APINO Shopper | | 18,20 |
| HERKULES padlock | | 10,85 |
| Passive light | | ● |
| Frame colour | | |
| Anthracite | | ● |

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge