



Company:

Basic price: **S-Eco 2 without drum brake** € 499,-

Basic price: **S-Eco 2 with drum brake** € 698,-

Customer:

Max. capacity 125 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 370 400 430 460 490 520	<input type="radio"/> 38,-
Seat height 500 mm	<input checked="" type="radio"/>
Seat depth 420 mm	<input checked="" type="radio"/>
Seat cushion Seat width 370 - 490 mm	24,-
Seat cushion Seat cushion 520 mm	30,-
Seat upholstery, Nylon	<input checked="" type="radio"/>
Sidepanels	
Sidepanel removable with arm support, swing-away, elevating	<input checked="" type="radio"/>
Sidepanel elevating left right Piece	108,20
Sidepanel elevating with foldable arm pads left right Piece	123,50
Leg support	
Swing-away leg support with individual footplates	<input checked="" type="radio"/>
Lower leg support elevating links rechts Piece	107,60
Amputee support links rechts Piece	157,-
Footplate with heel strap	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8 "	<input checked="" type="radio"/>
Transit version, Rear wheels 12,5" PU	210,-
Front wheels 7" x 30 mm	<input checked="" type="radio"/>
PU-tyres solid rubber	<input checked="" type="radio"/>
Pneumatic tyres	79,-
Spoke protection, transparent 24" left right Piece	145,-
Wheel back Pair	52,50
One arm drive only without drum brake left right	325,-
Hand rim rubber-coated Pair	110,-

▼ Please mark!

Back	€
Back support height 400 mm	<input checked="" type="radio"/>
Handgrip elevating, removable Pair **	152,40
Head support incl. fixation **	188,-
Back support upholstery, Nylon	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension Piece	23,-
Wheel lock for one-arm drive mechanical left right	88,60
Accessories	
Togo power push aid incl. mounting set	2.446,60
Mounting aid for adaptation	189,50
Anti-tip device Pair	87,50
Tilting aid	<input checked="" type="radio"/>
Safety belt	26,80
Crutch holder	23,70
Therapy tray	215,-
Oxygen/infusion holder and bracket	169,-
APINO Citybag	60,42
APINO Shopper	18,20
Nachrüstsatz Transit	210,-
Passive Beleuchtung	<input checked="" type="radio"/>
Frame colour	
Anthracite metallic	<input checked="" type="radio"/>

** Head support and elevating handgrips cannot be mounted together

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge