



Company:

Customer:

Date:

Basic price: **Neo XXL**
Max. capacity 200 kg

€ 8.090,-

▼ Please mark!

Seat	€
Seat width 450 mm effective 400 - 510 mm	●
Seat width adjustable 520 - 600 mm	205,-
Seat height 500 mm	●
Seat depth adjustable 420 - 510 mm	●
Seat cushion	●
Seat angle adjustable (-2° bis 19°) electrically mandatory with additional option R-Net control	1.390,-
Seat angle adjustable (-2° bis 19°) mechanically	●
Standard seat black	●
Comfort seat mandatory with extended seating frame	545,-
Extended seating frame	225,-
Recaro seat unit seat and back support not electrically adjustable	1.850,-
Back	
Back support height 440 mm	●
Back support angle adjustment (-5° bis 75°) electrically mandatory with additional option R-Net control	1.590,-
Back support angle adjustment (-5° bis 75°) mechanically	●
Back support adjustable by strap system	290,-
Head support incl. fixation	285,-
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Leg support adjustable mechanically	●
Leg support elevating with length adjustment, mechanically left right mandatory with calf supportn piece	405,-
Calf support for horizontally elevating leg support piece	120,-
Electrically elevating leg support mandatory with additional option R-Net control left right piece	990,-
Footplate angle adjustable	●
Motor and Battery	
Max. range 40 km	●
Max. Speed 6 km/h	●
Max. Speed 10 km/h incl. light-kit	595,-

▼ Please mark!

Battery 2x 12V / 60 Ah	●
Charger 24V / 8 Ah	●
Lighting	€
Front and rear lights, flashing lights	●
Control device	
VR-2 control	●
R-Net Can Bus mandatory for electrical adjustments	860,-
Control device, mounted left right	○
VR-2 control device for assistant complete with fixation	508,-
R-Net control device for assistant complete with fixation	718,-
Software kit: Software, cable (E)	490,-
Tyres	
Front and rear wheel single suspension	●
Rear wheels 3.00 - 8 pneumatic tyres	●
Front wheels 3.00 - 4 pneumatic tyres	●
Airless tyres RolkoFill	299,-
Rear mudguard piece	145,-
Accessories	
Anti-tip device	●
Transport attachment	●
Safety belt	71,20
Crutch holder	56,-
Rear mirror only with bracket	49,90
Bracket for rear mirror	69,-
Rear marking plate	75,-
PETITE Universal bag	9,30
BICO gripping pliers	26,61
APINO Citybag	60,42
APINO Shopper	18,20
Tool set	●
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge