



Company:

Customer:

Date:

Basic price: **Taiga**
Max. capacity 170 kg

€ 7.590,-

▼ Please mark!

Seat	€
Seat width 415-615 mm	●
Seat height 500 mm	●
Seat depth adjustable 410 - 560 mm	●
Seat cushion adjustable in width	●
Seat angle adjustable (-2° bis +10°) mechanically	●
Seat angle adjustable (-2° bis +20°) electrically	1.125,70
Standard seat black	●
Comfort seat 390 mm 440 mm 490 mm	530,50
Back	
Back support height 530 mm	●
Back support comfort seat 630 - 670 mm	●
Back support angle adjustable from 90° to 120°, mechanically	●
Back support angle adjustable from 90° to 135°, electrically	1.225,70
Back support adjustable by strap system not with comfort seat	●
Head support incl. fixation	144,-
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elevating, mechanically left right unit price	107,60
Leg support elevating, electrically left right unit price	458,35
Footplate angle adjustable	●
Motor and Battery	
Max. range approx. 35 km	●
Speed 6 km/h	●
Batteries 2x 12V/75 Ah	●
Battery charger 24V/8 Ah Standard	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

▼ Please mark!

Lighting	€
Front and rear lights, flashing lights	●
Control device	
swing-away control device, adjustable in length	●
VR-2 up to 2 electrical adjustments	●
VR-2 control device for assistant complete with fixation	508,-
R-Net up to 4 electrical adjustments swing-away control device, adjustable in length	2.490,-
R-Net control device for assistant complete with fixation	718,-
Control device mounted left right	○
Software kit: software, cable (E)	504,70
Tyres	
Rear wheel suspension	●
Rear wheel 14"	●
Front wheel 10"	●
Anti-tip device	●
Mudguard at rear wheels	●
Accessories	
Safety belt	30,90
Rear mirror only with bracket	51,40
Bracket for rear mirror	71,10
Kerb climber	265,70
Crutch holder	56,-
APINO Citybag	60,42
APINO Shopper	18,20
Tool set	●
Frame color	
Silver	●

● Basic equipment ○ Optional, no extra charge