



Company:

Basic price: **Terra**

€ 4.620,-

Customer:

Max. capacity 136 kg

Date:

▼ Please mark!

Seat	€
Seat width adjustable 410 - 500 mm	●
Seat height 540 mm	●
Seat depth adjustable 400 - 450* mm *standard	●
Seat cushion	●
Seat angle adjustable (-2° bis +10°) mechanically	●
Standard seat black	●
Comfort seat	530,50
Back	
Back support height 550 mm	●
Back support angle adjustable from 90° to 120°, mechanically	●
Back support adjustable by strap system	●
Head support incl. fixation	144,-
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elevating, mechanically left right unit price	107,60
Footplate angle adjustable	●
Width backdrop, adjustable	113,30
Motor and Battery	
Max. range approx. 32 km	●
Speed 6 km/h	●
Batteries 2x 12V / 55Ah	●
Battery charger 24V / 5Ah	●
Lighting	
Front and rear lights, flashing lights	●

▼ Please mark!

Control device	€
swing-away control device, adjustable in length	●
Control device mounted left right	○
VR-2 up to 2 electrical adjustments	●
VR-2 control device for assistant complete with fixation	508,-
Software kit: software, cable (E)	504,70
Tyres	
Rear wheel suspension	●
Rear wheel 3.00 - 8. PU	●
Front wheel 230 x 70 mm, PU	●
Mounting PU to pneumatic tires	79,-
Accessories	
Anti-tip device	●
Transport attachment	●
Kerb climber	265,70
Safety belt	●
Crutch holder	56,-
Rear mirror only with bracket	51,40
Bracket for rear mirror	71,10
APINO Citybag	60,42
APINO Shopper	18,20
Tool set	●
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge