



Company:

Customer:

Date:

Basic price: **Triplex**
Max. capacity 136 kg

€ 4.820,-

▼ Please mark!

Seat	€
Seat width adjustable 410 - 500 mm	●
Seat height 540 mm	●
Seat depth adjustable 400 - 450 mm	●
Seat cushion	●
Seat angle adjustable (-2° bis +10°) mechanically	●
Seat angle adjustable (-2° bis +20°) electrically	1.225,70
Standard seat Black	●
Comfort seat	530,50
Back	
Back support height 550 mm	●
Back support angle adjustable 90° to 120° mechanically	●
Back support infinitely adjustable 90° to 135° electrically	1.225,70
Back support adjustable by strap system	●
Head support incl. fixation	144,-
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elevating, mechanically left right unit price	107,60
Leg support, electrically left right unit price	458,35
Footplate angle adjustable	●
Width backdrop, adjustable	113,30
Motor and Battery	
Max. range approx. 32 km	●
Speed 6 km/h	●
Batteries 2x 12V / 55 Ah	●
Standard charger 24V / 5Ah	●

▼ Please mark

Lighting	€
Front and rear lights, flashing lights	●
Control device	
swing-away control device, adjustable in length	●
Control device, mounted left right	○
VR-2 up to 2 electrical adjustments swing-away control device, adjustable in length	●
VR-2 control device for assistant complete with fixation	508,-
R-Net up to 4 electrical adjustments swing-away control device, adjustable in length	2.490,-
R-Net control device for assistant complete with fixation	718,-
Software kit: Software, cable (E)	504,70
Tyres	
Rear wheel suspension	●
Rear wheel 3.00 - 8, PU	●
Front wheel 2.80/2.50 - 4, PU	●
Mounting PU to pneumatic tyres	79,-
Accessories	
Anti-tip device	●
Transport attachment	●
Kerb climber	265,70
Safety belt	●
Crutch holder	56,-
Rear mirror only with bracket	51,40
Bracket for rear mirror	71,10
APINO Citybag	60,42
APINO Shopper	18,20
Tool set	●
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge