



ORDER FORM

ECO XXL

Standard wheelchair

Company:

Customer:

Date:

Basic price: **Seat width 600 mm with hub brake** € 1.980,-

Basic price: **Seat width 650 mm with hub brake** € 1.980,-

Basic price: **Seat width 700 mm with hub brake** € 1.980,-

Max. capacity 250 kg

▼ Please mark!

Seat	€
Seat width in mm 600 650 700	<input type="radio"/>
Seat height 500 - 550 mm	<input checked="" type="radio"/>
Seat depth 480 mm	<input checked="" type="radio"/>
Seat cushion strengthened	<input checked="" type="radio"/>
Back support	
Back support height 460 - 535 mm	<input checked="" type="radio"/>
Back support adjustable by strap system in mm 600 650 700	203.50
Handgrip elevating, removable	<input checked="" type="radio"/>
Sidepanels	
Sidepanels swing away	<input checked="" type="radio"/>
Sidepanel elevating left right piece	129.70
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Footplate elevating	<input checked="" type="radio"/>
Footplate with heel strap	<input checked="" type="radio"/>
Lower leg support elevating left right	135.50
Amputee support left right	172.70

▼ Please mark!

Tyres	€
Rear wheels 24" PU strengthened axis	<input checked="" type="radio"/>
Front wheels 7" PU	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Drum brake	<input checked="" type="radio"/>
Accessories	
Tilting aid	<input checked="" type="radio"/>
Anti-tip device Pair	96.30
Crutch holder	<input checked="" type="radio"/>
Stabilization tube	<input checked="" type="radio"/>
Back support with bag	<input checked="" type="radio"/>
Calf strap	<input checked="" type="radio"/>
Oxygen/infusion holder and bracket	185.90
APINO Shopper	35.50
Passive light	<input checked="" type="radio"/>
Frame color	
Anthracite	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge