



ORDER FORM

PYRO START PLUS SL

Light weight wheelchair

Company:

Basic price: **Pyro Start Plus SL without drum brake** € 1.558.-

Basic price: **Pyro Start Plus SL with drum brake** € 1.817.-

Customer:

Max. capacity 150 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 430 460 490 520	○
Seat height front 525 mm adjustable from 525 - 600 mm in 25 mm steps	●
Seat depth 460 mm adjustable from 460 - 520 mm in 20mm steps.	●
Seat cushion Seat width 460 - 480 mm	33,-
Seat cushion Seat width 500 - 520 mm	37,40
Seat upholstery, Nylon	●
Back support	
Back support with lumbar curve	●
Back support height 460 mm adjustable 460 - 500 mm in 20 mm steps	●
Back support adjustable by strap-system	●
Handgrip elevating*, removable pair	167,70
Head support* incl. fixation	206,80
Back support upholstery, Nylon	●
Sidepanels	
Sidepanel removable with arm support, swing-away, elevating	●
Sidepanel elevating left right piece	119,10
Sidepanel elevating with foldable arm pads left right piece	135,90
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elongation 100 mm pair	67,-
Lower leg support elevating left right piece	118,40
Amputee support left right piece	172,70
Footplates with heel straps	●

* Head support and handgrips elevating cannot be mounted together

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	●
Transit version, rear wheels 12" PU	231,-
Front wheels 8" x 1 1/4", PU	●
Front wheels 200 x 50 mm PU pair	67,-
PU-tyres	●
Pneumatic tyres	86,90
Wheel base extension	●
Adapter block	●
Spoke protection, transparent 24" pair	159,50
One arm drive only without drum brake left right	357,50
Hand rim rubber-coated 24" pair	173,80
Wheel locks	
Knee lever wheel lock push to lock	●
Wheel lock lever extension piece	25,30
Knee lever wheel lock for one-arm drive left right	97,50
Accessories	
Togo power push aid incl. mounting set	2.676.-
Mounting aid for Togo adaptation	208,50
Retrofit kit Transit	231.-
Anti-tip device pair	96,30
Tilting aid	●
Safety belt	29,50
Crutch holder	26,10
Therapy tray	236,50
Oxygen/infusion holder and bracket	185,90
APINO Shopper	35,50
HERKULES padlock	25,50
Passive light	●
Frame colour	
Anthracite	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge