



ORDER FORM

S-VR Multifunctional wheelchair

Company:

Basic price: **S - VR without drum brake**

€ 1.668,-

Customer:

Basic price: **S - VR with drum brake**

€ 1.927,-

Date:

Max. capacity 120 kg

▼ Please mark!

Seat	€
Seat width in mm 400 440 480	<input type="radio"/>
Seat height 480 mm	<input checked="" type="radio"/>
Seat depth 420 mm	<input checked="" type="radio"/>
Seat angle adjustable 4°	<input checked="" type="radio"/>
Seat cushion 50 mm height	from 26,40
Back	
Back support height 430 / 850 mm	<input checked="" type="radio"/>
Back support angle -4° - +30° / +82°* adjusted mechanically	<input checked="" type="radio"/>
Neck holder adjustable in height, upholstered, removable	<input checked="" type="radio"/>
Back support upholstery, vinyl black	<input checked="" type="radio"/>
Sidepanels	
Sidepanels removable	<input checked="" type="radio"/>
Arm support elongated	<input checked="" type="radio"/>

▼ Please mark!

Leg support	€
Swing-away leg support with individual footplates plastic, removable, horizontally adjustable	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8" PU	<input checked="" type="radio"/>
Front wheels 200 x 50 mm PU	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Accessories	
Anti-tip device pair	<input checked="" type="radio"/>
Safety belt	29,50
Therapy tray	236,50
Fold blocking anti-theft device	<input checked="" type="radio"/>
Frame colour	
Anthracite	<input checked="" type="radio"/>

* Conversion not at factory

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge