



**Company:**

Basic price: **Terra**  
Max. capacity 136 kg

€ 5.313,-

**Customer:**

**Date:**

▼ Please mark!

| <b>Seat</b>   | €      |
|---|--------|
| Seat width adjustable 410 - 500 mm                                      | ●      |
| Seat height 540 mm  | ●      |
| Seat depth adjustable 400 - 450 mm                                      | ●      |
| Seat depth basic setting 450 mm   | ●      |
| Seat cushion  | ●      |
| Seat angle adjustable (-2° bis +10°)   mechanically                     | ●      |
| Standard seat   black   | ●      |
| Comfort seat  | 583,60 |
| <b>Back</b>   |        |
| Back support height 550 mm  | ●      |
| Back support angle adjustable in 4 steps from 90° to 120°, mechanically | ●      |
| Back support adjustable by strap system                                 | ●      |
| Head support   incl. fixation   | 158,40 |
| <b>Sidepanels</b>   |        |
| Adjustable in height and lateral  | ●      |
| <b>Leg support</b>  |        |
| Swing-away leg support with anatomically shaped handle                  | ●      |
| Leg support horizontal and width adjustable<br>left right   unit price  | 118,40 |
| Footplate angle adjustable  | ●      |
| Width backdrop, adjustable  | 124,70 |
| <b>Motor and Battery</b>  |        |
| Max. range approx. 32 km  | ●      |
| Speed 6 km/h  | ●      |
| Batteries 2 x 12V / 55Ah  | ●      |
| Battery charger 24V / 5A  | ●      |

▼ Please mark!

| <b>Lighting</b>  | €      |
|--|--------|
| Light system incl. indicators                              | ●      |
| <b>Control device</b>                                      |        |
| Swing-away control device, adjustable in length            | ●      |
| Control device mounted left right                          | ○      |
| VR-2 control device for assistant   complete with fixation | 995,50 |
| Software kit: software A + cable                           | 572,-  |
| <b>Tyres</b>   |        |
| Rear wheel suspension                                      | ●      |
| Rear wheel   3.00 - 8. PU                                  | ●      |
| Front wheel   230 x 70 mm, PU                              | ●      |
| Mounting PU to pneumatic tyres                             | 181,50 |
| <b>Accessories</b>   |        |
| Anti-tip device  | ●      |
| Transport attachment sideways                              | ●      |
| Kerb climber   | 292,30 |
| Safety belt  | ●      |
| Crutch holder  | 61,60  |
| Rear mirror   mountable only with bracket                  | 56,60  |
| Bracket for rear mirror                                    | 78,30  |
| APINO Shopper  | 35,50  |
| Tool set   | ●      |
| <b>Frame color</b>   |        |
| Silver   | ●      |

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment    ○ Optional, no extra charge