



ORDER FORM

CENTURO S2 4-Wheel Scooter

Company:

Basic price: **Centuro S2** | 6 km/h

€ 2.793.-

Basic price: **Centuro S2** | 10 km/h

€ 3.373.-

Customer:

Max. payload 160 kg

Date:

▼ Please Mark !

Technical Data	€
Range approx. 35 km	●
Ground clearance 80 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21%	●
Reversing width 1630 mm	●
Tyres 10" pneumatic	●
Sprung suspension front and rear	●
Brakes electromagnetic	●
Mechanical parking brake	●
Handlebars ergonomic	●
Adjustable steering column	●
Motor and electrical equipment	
Motor 6 km/h, 550 W	●
Motor 10 km/h, 600 W	●
Battery 2 x 12 V/55 Ah	●
Battery charger 5 A	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Lighting according to German road traffic regulations	●
Authorized for public transportation	●
Electronic speed reduction when cornering (only for 10 km/h)	●

▼ Please Mark !

Seat and back support	€
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●
Accessories	
Shopping basket removable	●
Floor mat	●
Rear mirror right	●
Rear mirror holder right	●
Rear mirror left*	19,80
Rear mirror holder left	2,90
Crutch holder for 2 Crutches Mounting at the seat frame	106,90
License plate holder required for >6 km/h **	26,60
Rain cover	101,70
Frame color	
White Anthracite	○

*just mountable with rear mirror holder

** In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge