



ORDER FORM

CENTURO S7 | S7 MAX | S7 XXL

4-Wheel Scooter

Company:

Customer:

Date:

Basic price: **Centuro S7** € 4.533,-
 Max. payload 160 kg

Basic price: **Centuro S7 MAX** € 4.663,-
 Max. payload 160 kg

Basic price: **Centuro S7 XXL** € 5.046,-
 Max. payload 250 kg

▼ Please Mark !

	S7	S7 MAX	S7 XXL
Technical Data			
Motor 15 km /h	●	X	X
Motor 10 km /h	X	●	●
Range approx. 45 km	●	●	X
Range approx. 35 km	X	X	●
Battery 2 x 12 V / 80 Ah		●	
Battery charger 8 A		●	
Tyres 13" pneumatic		●	
Ground clearance 130 mm		●	
Obstacle height 100 mm		●	
Climbing ability 12° 21 %		●	
Reversing width 2300 mm		●	
Standard equipment			
Sprung suspension front and rear		●	
Brakes electromagnetic		●	
Mechanical parking brake		●	
Handlebars ergonomic		●	
Adjustable steering column		●	
USB-Port		●	
Cockpit with LCD-Display		●	
Horn		●	
Rear mirror right and left		●	
Lighting according to German road traffic regulations		●	
Electronic speed reduction when cornering		●	

▼ Please Mark !

	S7	S7 MAX	S7 XXL
Seat and back support			
Swing-away seat, height-adjustable			
- Back support foldable and angle adjustable			
- Arm support foldable and angle adjustable	●	●	X
- Seat rails for adjusting the distance to the steering column			
XL Swing-away seat, height-adjustable	€	€	●
- Back support foldable and angle adjustable	165,-	165,-	
- Arm support foldable and angle adjustable	*	*	
- Seat rails for adjusting the distance to the steering column			
Accessories			
Shopping basket detachable		●	
Floor mat		●	
Cup holder		●	
Crutch holder for 1 Crutch	€	€	
Mounting on the multifunction holder	135,10		
Crutch holder for 2 Crutches	€	€	
Mounting at the seat frame	106,90		
Holder for cross-foldable rollators	€	€	
Mounting on the multifunction holder	237,60		
Holder for oxygen bottle	€	€	
Mounting on the multifunction holder	51,50		
Holder for transport box	€	€	
Mounting on the multifunction holder	53,50		
Transport box with lock	€	€	
	218,90		
License plate holder required for >6 km/h **	€	€	
	26,60		
Rain cover	€	€	
	101,70		
Frame color			
White	Anthracite		○

* Surcharge. Option price less price of standard equipment.
 This price differs from the unit price.

** In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge X not an Option