



# ORDER FORM

## PROTEGO SL Multifunctional wheelchair

**Company:**

Basic price: **Protego SL without drum brake** € 2.981,-

Basic price: **Protego SL with drum brake** € 3.178,-

**Customer:**

Max. capacity 150 kg

**Date:**

▼ Please mark!

Seat	€
Seat width in mm <b>440</b> - 465 <b>490</b> - 515	<input type="radio"/>
Seat height 520 - 600 mm   with seat cushion	<input checked="" type="radio"/>
Seat depth 480 - 520 mm	<input checked="" type="radio"/>
Seat angle adjustable 0 - 15°	<input checked="" type="radio"/>
Medial support	92,70
Seat upholstery Dartex, night grey	<input checked="" type="radio"/>
<b>Back</b>	
Back support   ergonomic shape with lateral support	<input checked="" type="radio"/>
Back support height 540 - 650 mm   with seat cushion	<input checked="" type="radio"/>
Back support angle 0 - 44°	<input checked="" type="radio"/>
Trunk support   swing away    left    right	94,80
Head support   incl. fixation	<input checked="" type="radio"/>
Back support upholstery Dartex, night grey	<input checked="" type="radio"/>
Handgrip elevating, removable	<input checked="" type="radio"/>
<b>Sidepanels</b>	
Sidepanels adjustable in height and depth, removable, with Dartex arm support	<input checked="" type="radio"/>
Hemiplegic arm support    left    right	298,70
<b>Leg support</b>	
Leg support elevating and angle adjustable with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU   incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" x 45 mm PU	<input checked="" type="radio"/>
Spoke protection, transparent 24"	145,-
<b>Wheel locks</b>	
Knee lever wheel lock   push to lock	<input checked="" type="radio"/>
Wheel lock lever extension	19,90
<b>Accessories</b>	
Togo power push aid incl. mounting set	2.523,-
Mounting aid for adaptation	189,50
Tilting aid	<input checked="" type="radio"/>
Safety belt	26,80
Therapy tray	215,-
Oxygen/infusion holder*    110 mm    160 mm	169,-
APINO Citybag	60,42
APINO Shopper	18,20
Passive light	<input checked="" type="radio"/>
<b>Frame colour</b>	
Anthracite	<input checked="" type="radio"/>

\* Cannot be mounted together with Togo

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment     Optional, no extra charge