



ORDER FORM

CENTURO S4 4-WHEEL SCOOTER

Company:

Basic price: **Centuro S4** | 6 km/h

€ 3.568.-

Customer:

Basic price: **Centuro S4** | 10 km/h

€ 4.311.-

Date:

Max. payload 160 kg

▼ Please Mark !

Technical Data	€
Range approx. 35 km	●
Ground clearance 115 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21%	●
Reversing width 1690 mm	●
Tyres 12" pneumatic	●
Sprung suspension front and rear	●
Brakes electromagnetic	●
Mechanical parking brake	●
Handlebars ergonomic	●
Adjustable steering column	●
Motor and electrical equipment	
Motor 6 km/h, 550 W	●
Motor 10 km/h, 700 W	●
Battery 2 x 12 V/55 Ah	●
Battery 2 x 12 V/80 Ah	288,30
Battery charger 5 A	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Lighting according to German road traffic regulations	●
Electronic speed reduction when cornering (only for 10 km/h)	●

▼ Please Mark !

Seat and back support	€
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●
XL Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	165,-
Accessories	
Shopping basket removable	●
Floormat	●
Rear mirror right	●
Rear mirror holder right	●
Rear mirror left*	20,50
Rear mirror holder left	3,-
Crutch holder for 1 Crutch Mounting on the multifunction holder	140,50
Crutch holder for 2 Crutches Mounting at the seat frame	111,10
Holder for cross-foldable rollators Mounting on the multifunction holder	247,10
Holder for oxygen bottle Mounting on the multifunction holder	53,10
Holder for transport box Mounting on the multifunction holder	55,64
Transport box with lock	227,60
License plate holder required for > 6 km/h **	27,60
Rain cover	105,70
Frame color	
White Anthracite	○

*just mountable with rear mirror holder

** In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge